

ACUVUE®

# YOUR QUICK AND EASY GUIDE TO FIT SUCCESS



Product images for illustrative purposes only



**Unique  
PUPIL OPTIMISED  
DESIGN<sup>1</sup>**

**ACUVUE®  
MULTIFOCAL  
Fit guide**

**Fit success  
& patient  
satisfaction**

## DESIGNED FOR SUPERIOR VISUAL PERFORMANCE.\*<sup>1</sup>

Now available as both Daily Disposable  
and Reusable contact lenses.


ACUVUE®  
**MULTIFOCAL**  
WITH PUPIL OPTIMISED DESIGN

\* Compared to prior JJV multifocal design; technology optimised for both the parameters of refractive error and add power for a multitude of viewing distances and light levels



The only brand with 100% of parameters  
**optimised by age & refraction\*\*1**

## PUPIL OPTIMISED DESIGN



For illustrative purposes only. Pupil area can vary by ~20% at a given luminance.<sup>3</sup>

## IN-BUILT PRECISION

ACUVUE® MULTIFOCAL PORTFOLIO with PUPIL OPTIMISED DESIGN provides a more **PRECISE FIT**: Hybrid Back Curve Technology **better matches the shape of the natural eye** to help keep the lens' optical design in the **right place**.<sup>1</sup>



Product images for illustrative purposes only

**Every parameter is designed to match different pupil sizes and provide the best balance of vision for that age and refraction<sup>1</sup>**

\*\* Compared to competitor's designs; technology optimised for both the parameters of refractive error and add power.  
 † Across the power range of +6.00D to -9.00D.



# Use the **Fit Guide** for a more precise fit<sup>1</sup>

## INITIAL LENS SELECTION

**1 Determine the Best Vision Sphere (BVS)**  
 In the trial frame, confirm the least minus spherical prescription that provides the best distance VA<sup>†</sup>

**2 Determine the sensory dominant eye**  
 The +1.00D blur test recommended with the BVS in the trial frame rather than sighting methods.

**3 Determine the lowest ADD based upon the patient's needs**  
 With the BVS in the trial frame, now determine the lowest ADD required to achieve good near vision.  
**Top Tip:** Start with 0.50D less than the spectacle ADD and if necessary, increase in 0.25D steps until required near vision is achieved.

**4 Select lens based on following tables**

Spectacle ADD	Initial Lens Selection		Enhance Distance		Enhance Near	
	Dominant Eye	Non-Dominant Eye	Dominant Eye	Non-Dominant Eye	Dominant Eye	Non-Dominant Eye
+0.75 to +1.25			Use a spherical ACUVUE <sup>®</sup> lens			 & give extra +0.25D to dist. Rx
+1.50 to +1.75						 & give extra +0.25D to dist. Rx
+2.00 to +2.50				 & give extra +0.25D to dist. Rx		 & give extra +0.25D to dist. Rx

<sup>†</sup> Proceed if astigmatism is less than or equal to 0.75DC. <sup>9</sup> Apply vertex distance correction if greater than +/- 4.00D.



# Increase patient and practice success by following the Fit Guide

## ACUVUE® MULTIFOCAL PORTFOLIO WITH PUPIL OPTIMISED DESIGN TECHNOLOGY:



Offers a more **personalised solution** for your patients.\*\*1



Provides a more **precise fit** to help keep your patient's optics in the right place and the right shape<sup>1</sup>

LENS DETAILS	1-DAY ACUVUE® MOIST MULTIFOCAL	ACUVUE® OASYS MULTIFOCAL 2-WEEKLY	ACUVUE® OASYS MAX 1-Day MULTIFOCAL
Material	etafilcon A	senofilcon A	senofilcon A
Diameter	14.3 mm	14.3 mm	14.3 mm
Base curve	8.4 mm	8.4 mm	8.4 mm
Technology	Embedded PVP*/LACREON® Technology	Embedded PVP*/HYDRACLEAR® PLUS Technology	TearStable™ Technology OptiBlue™ Light Filter‡
UV blocker*	Class 2	Class 1	Class 1
Dk/t*	25.5 x 10 <sup>-9</sup>	147 x 10 <sup>-9</sup>	147 x 10 <sup>-9</sup>
Visibility tint	Yes	Yes	Yes   Blue-green§2
Sphere	-9.00D to +6.00D (0.25D steps)	-9.00D to +6.00D (0.25D steps)	-9.00D to +6.00D (0.25D steps)
ADD	LOW +0.75D to +1.25D MID +1.50D to +1.75D HIGH +2.00D to +2.50D	LOW +0.75D to +1.25D MID +1.50D to +1.75D HIGH +2.00D to +2.50D	LOW: +0.75D to +1.25D MID: +1.50D to +1.75D HIGH: +2.00D to +2.50D

# Oxygen transmissibility at centre of a -3.00D lens using boundary-corrected, edge-corrected Dk values. Units: (cm/sec) (ml O2/ml x mm Hg) at 35°C. Dk determined via polarographic method.

+PVP=polyvinylpyrrolidone.



Visit the **ACUVUE® Multifocal Fitting Calculator** for quick & easy contact lens fitting & lens selection



\*\* Compared to competitor's designs; technology optimised for both the parameters of refractive error and add power. \* All ACUVUE® contact lenses have Class 1 or Class 2 UV-blocking to help provide protection against transmission of harmful UV radiation to the cornea and into the eye. UV-absorbing contact lenses are NOT substitutes for protective UV absorbing eyewear such as UV-absorbing goggles or sunglasses because they do not completely cover the eye and surrounding area. UV transmission measured with -1.00D lens. † Filtering of High Energy Visible (HEV) light by contact lenses has not been demonstrated to confer any systemic and/or ocular health benefit to the user. The Eye Care Professional should be consulted for more information. ‡ ACUVUE® OASYS MAX 1-Day has a unique blue-green appearance as a result of the combination of the blue-violet/high energy visible [HEV] light filter and the blue handling tint.

1. JJV Data on File 2022 ACUVUE® PUPIL OPTIMISED DESIGN TECHNOLOGY: JVCV Contact Lenses, Design Features, and Associated Benefits.
2. JJV Data on File 2022. TearStable™ Technology Definition.

For more information on proper wear, care and safety, please consult the Instructions for Use or visit our J&J website [www.jnjvisionpro.co.uk](http://www.jnjvisionpro.co.uk).

© Johnson & Johnson Medical Ltd 2022. ACUVUE®, ACUVUE® MOIST, ACUVUE® OASYS, ACUVUE® OASYS MAX 1-Day MULTIFOCAL, LACREON®, HYDRACLEAR®, TearStable™, OptiBlue™ are registered trademarks of Johnson & Johnson. PP2022MLT7062

**ACUVUE®**  
**MULTIFOCAL**  
WITH PUPIL OPTIMISED DESIGN