



# Increase patient and practice success by following the Fit Guide

## ACUVUE® MULTIFOCAL PORTFOLIO with PUPIL OPTIMISED DESIGN TECHNOLOGY:



Offers a more **personalised solution** for your patients.\*\*<sup>3</sup>



Delivers **consistent success** when fitted to a wide range of presbyopes, **regardless of vision correction needs.**<sup>3</sup>

LENS DETAILS <sup>4</sup>	1-DAY ACUVUE® MOIST MULTIFOCAL	ACUVUE® OASYS MULTIFOCAL 2-WEEKLY
Material	etafilcon A	senofilcon A
Diameter	14.3 mm	14.3 mm
Base curve	8.4 mm	8.4 mm
Wetting agent	Embedded PVP <sup>†</sup> /LACREON® Technology	Embedded PVP <sup>†</sup> /HYDRACLEAR® PLUS Technology
UV blocker <sup>‡</sup>	Class 2	Class 1
Dk/t <sup>§</sup>	25.5 x 10 <sup>-9</sup> §	147 x 10 <sup>-9</sup> §
Visibility tint	Yes	Yes
Sphere	-9.00D to +6.00D (0.25D steps)	-9.00D to +6.00D (0.25D steps)
ADD	LOW +0.75D to +1.25D MID +1.50D to +1.75D HIGH +2.00D to +2.50D	LOW +0.75D to +1.25D MID +1.50D to +1.75D HIGH +2.00D to +2.50D

§ Dk/t units: 10 (cm/sec) (mL O<sub>2</sub>/mL x mm Hg). All Dk values: Fatt units at 35°C, determined via polarographic method (boundary & edge corrected) (-3.00D lens).

† PVP=polyvinylpyrrolidone.



Visit the **ACUVUE® Multifocal Fitting Calculator** for quick & easy contact lens fitting & lens selection - [www.jnjvisioncare.co.uk/5minfit](http://www.jnjvisioncare.co.uk/5minfit)

‡ All ACUVUE® contact lenses have Class 1 or Class 2 UV-blocking to help provide protection against transmission of harmful UV radiation to the cornea and into the eye. UV-absorbing contact lenses are NOT substitutes for protective UV absorbing eyewear such as UV-absorbing goggles or sunglasses because they do not completely cover the eye and surrounding area. UV transmission measured with -1.00D lens.

# With four total lenses or less.

\*Compared to prior JJV multifocal design; technology optimised for both the parameters of refractive error and add power for a multitude of viewing distances and light levels.

\*\*Compared to competitor's designs; technology optimised for both the parameters of refractive error and add power.

1. JJV Data on File 2015. 1-DAY ACUVUE® MOIST MULTIFOCAL Fit and Performance.

2. JJV Data on File 2020. ACUVUE® OASYS MULTIFOCAL Fit and Performance Claims.

3. JJV Data on File 2021 ACUVUE® PUPIL OPTIMISED DESIGN TECHNOLOGY: JJVC Contact Lenses, Design Features, and Associated Benefits.

4. JJV Data on file 2015. 1-DAY ACUVUE® MOIST MULTIFOCAL lens details.

ACUVUE®, 1-DAY ACUVUE® MOIST, and LACREON® are trademarks of Johnson & Johnson. © Johnson & Johnson 2021. P2021MLT5140

ACUVUE®  
**MULTIFOCAL**  
WITH PUPIL OPTIMISED DESIGN

ACUVUE




# YOUR QUICK AND EASY GUIDE TO FIT SUCCESS

NEW



94% fit success<sup>#1</sup>




97% fit success<sup>#2</sup>

Product images for illustrative purposes only.



Unique  
PUPIL OPTIMISED DESIGN



ACUVUE® MULTIFOCAL  
Fit guide



Fit success & patient  
satisfaction


Designed for exceptional visual performance.\*<sup>3</sup>  
Now available as both Daily Disposable and Reusable contact lenses.

ACUVUE®  
**MULTIFOCAL**  
WITH PUPIL OPTIMISED DESIGN



The brand with 100% of parameters **optimised by age & refraction\*\*3**

### PUPIL OPTIMISED DESIGN



For illustrative purposes only. Pupil area can vary by ~20% at a given luminance.\*

### IN-BUILT PRECISION

ACUVUE® MULTIFOCAL PORTFOLIO with PUPIL OPTIMISED DESIGN provides a more PRECISE FIT: Hybrid Back Curve Technology **better matches the shape of the natural eye** to help keep your patient's optics in the **right place and the right shape**.<sup>3</sup>



**NEW**

**Same design, parameters and fit process across modalities so you can fit with success in a few easy steps<sup>1</sup>**



Use the **Fit Guide** for quick and easy fit success<sup>#1,2</sup>

### Initial Lens Selection

- Determine the Best Vision Sphere (BVS)**  
Perform a new subjective refraction<sup>o</sup> then remove cylindrical power and determine best distance VA with the least minus spherical refraction.<sup>oo</sup>
  - Determine the sensory dominant eye**  
+1.00D blur test recommended rather than sighting methods.
  - Determine the lowest ADD based upon the patient's needs.** With the distance BVS in place, determine the lowest functional reading ADD.
- Top Tip:** Over the spherical distance refraction, start with 0.50D less than the spectacle ADD and if necessary, increase in 0.25D steps until required near vision is achieved.

- Refer to the selection table based on the ADD, and determine initial trial lenses**

ADD	EYE	LENS SELECTION
+0.75D to +1.25D	Dominant Eye	LOW
	Non-dominant Eye	LOW
+1.50D to +1.75D	Dominant Eye	MID
	Non-dominant Eye	MID
+2.00D to +2.50D	Dominant Eye	MID
	Non-dominant Eye	HIGH

Allow for 10 minutes of real-world exposure (outside of the exam room) and proceed to trial, unless patient expectations or required standards for driving are not met.

**10 min**

- If an enhancement is needed, refer to the enhancement tables below - an over-refraction is not recommended**

ENHANCED DISTANCE VISION	ENHANCED NEAR VISION
SPHERICAL ACUVUE® lens	LOW
LOW	LOW +0.25D <sup>^</sup>
LOW	MID
MID	MID +0.25D <sup>^</sup>
MID	MID
MID +0.25D <sup>^</sup>	HIGH +0.25D <sup>^</sup>

<sup>^</sup> Add +0.25D to the distance power.

\*Across the power range of +6.00D to -9.00D.

† Up to 4 lenses total. <sup>o</sup> Proceed if astigmatism is less than 0.75DC. <sup>oo</sup> Apply vertex distance correction if greater than +/- 4.00D.